

Phonak Wireless Communication Portfolio



With today's technology, hearing aids have made speech understanding enjoyable, even in the most challenging listening situations. But sometimes, an extra boost can be useful. This is why Phonak offers a wide range of wireless accessories which boost your hearing aids performance and keep you connected. Whether on the phone, watching TV or in a noisy environment, we offer the optimal solutions to suit your everyday needs while benefitting from the exceptional performance of your hearing aids.

Control

Easy to use remote control for your hearing aids and an App for discreet access.



PilotOne II



RemoteControl App II¹

Phone

Talk on the phone at home and on your mobile phone with ease.



DECT II



EasyCall II

Connectivity

Wirelessly stream TV, music, and conversations in stereo quality.



ComPilot Air II



ComPilot II



TVLink II¹

Noise

Understand speech in noise and over distance in groups or one-on-one.



Roger Pen
Roger EasyPen



Roger Clip-On Mic



Roger receiver

Distance

Enhance one-to-one conversations over distance.



RemoteMic

WorkLife

Full participation in meetings



Roger Table Mic

¹ Requires Phonak ComPilot II or ComPilot Air II

Phonak has the right hearing solution for you – regardless of your type of hearing loss, lifestyle, personal preferences, age or budget.

Premium

Developed to satisfy those who desire and expect nothing but the best. Thanks to a complete set of trailblazing solutions that uniquely address all hearing challenges, these hearing aids support understanding in even the most difficult listening situations, for example at a busy social event.

Advanced

All the performance levels have the perfect balance between performance and price. An impressive offering of features that enhances hearing ability in a variety of listening environments, including those with moderate levels of background noise, for example at a restaurant.

Choose to enjoy better hearing that matches your lifestyle!

Next appointment:

Date _____

Time _____

This is your personal record of your solution. Please bring it with you to every appointment. Feel free to add notes of your own.

Standard

Ideally suited for people who expect value in state-of-the-art technology. A competitive feature set that has the capability to support the hearing demands of daily life with minimal to moderate levels of background noise, for example using the telephone in a public place.

Essential

The Essential products are designed for people looking for an affordable yet effective solution to enhance daily communication abilities, mostly in quiet or minimal background noise, for example one to one conversations.

Your hearing care professional:




Name _____

Date _____

¹ Requires Phonak ComPilot II or ComPilot Air II

Your audiogram



Hearing aid styles

Together with your hearing care professional, you can decide which hearing aid will best suit your needs and preferences.


ITE (In-the-Ear) ITE hearing aids are available as custom and extended wear models. Custom models are individually crafted to your ear canal to fit perfectly and comfortably in your ears, providing absolute discretion. Extended wear models are completely invisible, worn around the clock with no hassle.



RIC (Receiver-In-Canal) RIC models are small and slim, which makes them discreet. They literally disappear behind the ear and are very comfortable to wear.



BTE (Behind-the-Ear) BTE models are often recommended for people with more significant hearing loss or when dexterity is limited. They sit comfortably behind the ear, are easy to handle and very robust.



Available in a variety of hair, skin, fashion and fun colors.

Select your hearing goals



Premium (P)

Follow discussions in the car

Enjoy conversations at noisy social events

Advanced (A)

Hear voices from any direction

Enjoy fuller and richer music

Enjoy conversations at restaurants

Standard (S)

Converse with colleagues

Chat on the phone

Essential (E)

Enjoy personal discussions

Your recommended performance level _____

Notes _____

Phonak programs and features overview



What makes Phonak hearing aids so special is their unique operating system called AutoSense OS™. It senses your surroundings and adjusts your hearing aids every step of the way. AutoSense OS analyzes the sounds around you every 0.4 seconds and can identify whether you're in a noisy restaurant, car, concert hall or at home. It draws upon multiple features blending them to create over 200 distinct settings to precisely match your exact sound environment. Without having to adjust your hearing aids manually, you can enjoy unmatched hearing performance... wherever life takes you.

	P	A	S	E	
AutoSense OS™ programs	Comfort in Echo Reduces distracting echo, making listening in echoey environments such as large public buildings more comfortable.	•			
	Speech in Car Listening effort is reduced so you can enjoy the journey and the conversation.	•			
	Speech in Loud Noise Provides understanding of a single voice in front of you, for example while gathering with your family in a restaurant.	•			
	Music Offers a full and rich music experience so you can enjoy listening to your favorite music.	•	•		
	Comfort in Noise Uncomfortable loud sounds are reduced enhancing your listening comfort.	•	•	•	
	Speech in Noise Noise is reduced from the back allowing you to focus on the conversation in front of you.	•	•	•	•
	Calm Situation Creates a comfortable listening environment to simply relax or enjoy a one to one conversation.	•	•	•	•
Additional programs	Speech in Wind Allows you to enjoy conversations even on a windy day, for example while playing golf.	•			
	Speech in 360° Automatically selects the direction of speech, improving speech intelligibility without the need for you to face the speaker. Follow conversations for example in a business meeting around a large table or while driving.	•	•		
	Speech in Loud Noise Provides understanding of a single person speaking directly in front of you in a noisy crowd, for example while gathering with your family in a restaurant.	•	•		
Signature features	Acoustic phone Makes understanding on the phone easy by streaming the caller's voice to both ears, any time, any where. Talk on your home or mobile phone.	•	•	•	
	SoundRecover2 Makes high-pitched sounds easier to hear in all situations, for example birds singing, children's voices or the ring of a doorbell.	•	•	•	•
	WhistleBlock Eliminates annoying whistling from hearing aids so you can keep focused on the conversation without worry of interruption.	•	•	•	•
	NoiseBlock Reduces uncomfortable background noise such as traffic or the constant hum of a ceiling fan without affecting your ability to understand.	•	•	•	•

For a full list of Phonak program and feature benefits, please visit www.phonak.com.au or speak to your hearing care professional.