



ReSound GN

Quick reference guide

**ReSound LiNX<sup>2</sup>**™ and **ReSound ENZO<sup>2</sup>**™

The world's first full family of  
Smart Hearing™ aids.



DESIGN

## Custom Models

		Non-W	W	Options	LP: Mild-Moderate MP: Mild-Moderate	HP: Moderate - Severe UP: Severe - Profound	Matrices
	IIC	●		10A battery Faceplate color: anthracite, beige tinnitus sound generator, SoundShaper, Acceptance Manager			low power: 115/40
	CIC	●		10A battery VC, push button, tinnitus sound generator, SoundShaper, Acceptance Manager			
	MIH-S	●		10A battery VC, push button, tinnitus sound generator, SoundShaper, Acceptance Manager			
	MIH	●	●	10A, 312 or 13 battery Made for iPhone, Ear-to-Ear communication, push button, VC, telecoil, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager, 10A (non-wireless only)			low power: 115/40 medium power: 119/50 high power: 121/60 ultra power: 130/70
	ITC	●	●	312 battery Made for iPhone, directional microphones, telecoil, VC, push button, SoundShaper, Acceptance Manager, tinnitus sound generator			
	ITE HS	●	●	312 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			medium power: 119/50 high power: 121/60 ultra power: 130/70
	ITE FS	●	●	312 or 13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			

### Available Colors

ReSound custom models








ReSound LiNX<sup>2</sup>™ RIE and BTE



ReSound ENZO<sup>2</sup>™ BTE



## RIE/BTE Models

		Non-W	W	Options	LP: Mild-Moderate MP: Mild-Moderate	HP: Moderate - Severe UP: Severe - Profound	Matrices
	Mini RIE 61		●	312 battery <i>options:</i> Made for iPhone, Ear-to-Ear communication, directional microphones, programmable push button (PB or VC), tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">low power: 113/52 medium power: 116/58 high power: 122/65 ultra power: 129/75</p>
	RIE 62		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, multifunction button (VC, PB, & streaming), DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			
	Mini BTE 67		●	312 battery Made for iPhone, Ear-to-Ear communication, directional microphones, earhook or Thin Tube fitting, telecoil, push button, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">open: 122/52 closed: 120/56</p>
	Standard BTE 77		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, earhook or Thin Tube fitting, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">open: 127/51 closed: 124/57</p>
	Power BTE 88		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">power: 132/67</p>
	High Power BTE 88		●	13 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">power: 132/73</p>
	Super Power BTE 98		●	675 battery Made for iPhone, Ear-to-Ear communication, directional microphones, telecoil, VC, push button, DAI, tinnitus sound generator, SoundShaper, Synchronized Acceptance Manager			<p style="text-align: center;">power: 141/83</p>



## CONNECTIVITY

### ReSound wireless connections

Your patients can connect directly to all ReSound wireless accessories for direct audio streaming and discreet control over their listening experience - without using a neckloop. With a full ecosystem of 2.4 GHz wireless connectivity, your patients can hear everything better - from family members to the TV - with less effort.

Easily set up by pressing the pairing button on the accessory, open and close the battery doors of the hearing aids, and the accessory is automatically paired with wireless ReSound hearing aids and ready for wireless communication.

- **ReSound Micro Mic:** Significantly improves signal-to-noise ratio in noisy situations, delivers exceptional speech understanding and extends the hearing range to over 80 feet.
- **ReSound Multi Mic:** With the same hearing range of over 80 feet, doubles as a table microphone, connects with loop and FM systems\*\*, and has a mini-jack input for streaming audio from a computer or mobile device.
- **ReSound Unite Remote Control 2:** Adjust the volume, mute hearing aids, switch programs and see all settings at a glance.
- **ReSound Unite Phone Clip+:** Make crystal clear phone calls, mute background noise while talking, or stream music, podcasts or any other audio from a Bluetooth-enabled phone.
- **ReSound Unite TV 2:** Stream stereo sound from a TV, computer or music system directly to the hearing aids over a range of up to 20 feet. Includes a Dolby® adaptor for newer televisions.

### Airlink 2

With Airlink®2, the same wireless technology can also be used for programming. Imagine fitting without the use of any cables—at all!



\*\*FM receiver required

©2016 The GN ReSound Group, all rights reserved. Apple, iPad, iPhone and iPod touch are trademarks of Apple Inc, registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android and Google Play are registered trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.



ReSound Micro Mic



ReSound Multi Mic



ReSound Unite® Phone Clip+



ReSound Unite® Remote Control 2



ReSound Unite® TV Streamer 2

### Manufacturer according to FDA:

**ReSound North America**  
8001 Bloomington Freeway  
Bloomington, MN 55420  
1-800-248-4327  
resoundpro.com

### Manufacturer according to Health Canada:

**ReSound Canada**  
303 Supertest Road  
Toronto, Ontario M3J 2M4  
1-888-737-6863  
resoundpro.com



MK604322 Rev G 2016.07

# ReSound LiNX<sup>2</sup> technology levels

Functional Features	ReSound LiNX <sup>2</sup> ™ 9	ReSound LiNX <sup>2</sup> ™ 7	ReSound LiNX <sup>2</sup> ™ 5
Fully flexible programs	4	4	4
Ear-to-Ear communication	●	○	○
Synchronized push button	●	●	●
Synchronized volume control	●	●	●
SmartStart™	●	●	●
PhoneNow™	●	●	●
Comfort Phone™	●	●	●
Direct audio streaming (Made for iPhone)	●	●	●
ReSound Unite™ TV Streamer 2, Remote Control 2, Phone Clip+, Mini Microphone	●	●	●
ReSound Smart™ app	●	●	●
ReSound Control™ app (Phone Clip+ required)	●	●	●
<b>Audiological Features</b>			
WARP compression - number of channels	17	17	9
Environmental Classifier	●	●	●
Spatial Sense™	●		
Directionality™	●	○	○
Directional Mix Processor	●	●	●
-Adjustable directional mix	●	●	
Softswitching™	●	●	○
Adaptive Directionality™	●	○	○
Environmental Optimizer	●	○	
NoiseTracker™ II	●	○	○
Expansion	●	○	○
Windguard™	●	○	○
Sound Shaper	●	●	●
Low Frequency Boost (UP receiver only)	●	○	○
DFS Ultra™ II	●	●	●
-Music Mode™	●	●	●
Auto DFS™	●	●	●
Synchronized Acceptance Manager	●	●	●
Amplification strategy (WDRC/Semilinear/Linear - UP receiver only)	●	○	○
Tinnitus Sound Generator with Nature Sounds	●	●	●
<b>Fitting Features</b>			
Fitting software Aventa 3.8 or higher	●	●	●
Available gain handles	9	7	6
Onboard Analyzer™ II	●	●	●
In Situ Audiometry	●	●	●
Wireless fitting with Airlink™	●	●	●

○ Basic

○ Advanced

● Ultimate



	FULL PERFORMANCE			ALMOST INVISIBLE	SUPER POWER	PEDIATRICS	
	DESIGN PERFORMANCE	BTE PERFORMANCE	CUSTOM PERFORMANCE				
PREMIUM	 LiNX² 9 RIE 61	 LiNX² 9 RIE 62   BTE 77 / 88	 LiNX² 9 CIC   ITC   ITE	 LiNX² 9 IIC	 ENZO² 9 BTE 88 / 98	 Up Smart 9 BTE 77 / 88 / 98	
ADVANCED	 LiNX² 7 & 5 RIE 61	 LiNX² 7 & 5 RIE 62   BTE 77 / 88	 LiNX² 7 & 5 CIC   ITC   ITE	 Lex 8 Instant Fit   LiNX² 7 & 5 IIC	 ENZO² 7 & 5 BTE 88 / 98	 Up Smart 7 & 5 BTE 77 / 88 / 98	
ESSENTIAL		 Enya 4, 3 & 2 RIE 62   BTE 67 / 77 / 88	 Enya 4, 3 & 2 CIC   ITC   ITE	 Lex 4 Instant fit	 Magna 4 & 2 BTE 90		
CONNECTIVITY	 <b>Phone Clip+</b> Stream phone conversations and stereo sound from your Bluetooth enabled phone directly to your hearing instruments.		 <b>Multi Microphone</b> Significantly improves SNR in noisy situations, hearing range of 25 metres, doubles as a table microphone, connects with loop and FM systems and streams audio from a computer or music player.		 <b>TV Streamer 2</b> Stream sound from TV, computer, stereo or other audio devices directly to hearing instruments.		<b>Remote Control 2</b> Adjust volume monaurally and binaurally, mute hearing aids, change programs or switch between TV, radio, computer, etc.
APPS	 <b>ReSound Smart</b> Fine-tune hearing aid settings directly from your Smartphone. Create favourite programs, find lost hearing aids, control your Unite accessories, etc.		 <b>ReSound Control</b> Adjust volume and change programs, control your Unite accessories from your Smartphone, through the Phone Clip+.		 <b>ReSound Relief</b> Multisensory Tinnitus app that combines a variety of sounds, including customisable and dynamic soundscapes and guided relaxation exercises to help divert attention away from tinnitus.		



### Smart Hearing

A new kind of hearing aid that adapts to the way people want to live today, combining the best available sound quality with new possibilities to connect with their world. Smart Hearing aids are wireless.



### Wireless

ReSound wireless hearing aids stream sound using 2.4 GHz technology for:

- Crystal clear sound
- Stereo streaming
- No neckworn device
- Streaming at least 7 meters



### Tinnitus

Comprehensive tinnitus management solution to help relieve symptoms for those with tinnitus and hearing loss. Includes built-in tinnitus sound generator, nature sounds, apps and more



# ReSound Enya

## DEVICE SUMMARY

MODEL	STYLE	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY	
<b>Enya 4 Wireless</b>				
Enya 4 EY462-DRW	RIE	✓	312	
Enya 4 EY467-DW	Mini BTE	✓	13	
Enya 4 EY477-DW	Standard BTE	✓	13	
Enya 4 EY488-DW	Power BTE	✓	13	
Enya 4 EY4ITE-DW	ITE (MP)	✓	13	
	ITE (HP)	✓	13	
	ITE (UP)	✓	13	
Enya 4 EY4ITC-DW	ITC (LP)	✓	312	
	ITC (MP)	✓	312	
	ITC (HP)	✓	312	
	ITC (UP)	✓	312	
<b>Enya 3 Wireless</b>				
Enya 3 EY362-DRW	RIE	✓	312	
Enya 3 EY367-DW	Mini BTE	✓	13	
Enya 3 EY377-DW	Standard BTE	✓	13	
Enya 3 EY388-DW	Power BTE	✓	13	
Enya 3 EY3ITE-DW	ITE (MP)	✓	13	
	ITE (HP)	✓	13	
	ITE (UP)	✓	13	
Enya 3 EY3ITC-DW	ITC (LP)	✓	312	
	ITC (MP)	✓	312	
	ITC (HP)	✓	312	
	ITC (UP)	✓	312	
MODEL	STYLE	FTC	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY
<b>Enya 2 Wireless</b>				
Enya 2 EY267-DW	Mini BTE	✓	RC2 ONLY	312
Enya 2 EY277-DW	Standard BTE	✓	RC2 ONLY	13
Enya 2 EY288-DW	Power BTE	✓	RC2 ONLY	13
Enya 2 EY2ITE-DW	ITE (MP)	✓	RC2 ONLY	13
	ITE (HP)	✓	RC2 ONLY	13
	ITE (UP)	✓	RC2 ONLY	13
Enya 2 EY2ITC-DW	ITC (LP)	✓	RC2 ONLY	312
	ITC (MP)	✓	RC2 ONLY	312
	ITC (HP)	✓	RC2 ONLY	312

MODEL	STYLE	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY	
<b>Enya 4 Non Wireless</b>				
Enya 4 EY4ITE-D	ITE (MP)	-	13	
	ITE (HP)	-	13	
	ITE (UP)	-	13	
Enya 4 EY4ITC-D	ITC (LP)	-	312	
	ITC (MP)	-	312	
	ITC (HP)	-	312	
	ITC (UP)	-	312	
Enya 4 EY4CIC	CIC (LP)	-	10A	
	CIC (MP)	-	10A	
	CIC (HP)	-	10A	
	CIC (UP)	-	10A	
<b>Enya 3 Non Wireless</b>				
Enya 3 EY3ITE-D	ITE (MP)	-	13	
	ITE (HP)	-	13	
	ITE (UP)	-	13	
Enya 3 EY3ITC-D	ITC (LP)	-	312	
	ITC (MP)	-	312	
	ITC (HP)	-	312	
	ITC (UP)	-	312	
Enya 3 EY3CIC	CIC (LP)	-	10A	
	CIC (MP)	-	10A	
	CIC (HP)	-	10A	
	CIC (UP)	-	10A	
MODEL	STYLE	FTC	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY
<b>Enya 2 Non Wireless</b>				
Enya 2 EY2ITE-D	ITE (MP)	✓	-	13
	ITE (HP)	✓	-	13
	ITE (UP)	✓	-	13
Enya 2 EY2ITC-D	ITC (LP)	✓	-	312
	ITC (MP)	✓	-	312
	ITC (HP)	✓	-	312
	ITC (UP)	✓	-	312
Enya 2 EY2CIC	CIC (LP)	✓	-	10
	CIC (MP)	✓	-	10
	CIC (HP)	✓	-	10
	CIC (UP)	✓	-	10



# ReSound ENZO<sup>2</sup>

## DEVICE SUMMARY

MODEL	STYLE	MFI	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY
<b>ENZO<sup>2</sup> 9 Wireless</b>				
ENZO <sup>2</sup> 9 EN988	Power BTE	✓	✓	13
ENZO <sup>2</sup> 9 EN998	Super Power BTE	✓	✓	675
<b>ENZO<sup>2</sup> 7 Wireless</b>				
ENZO <sup>2</sup> 7 EN788	Power BTE	✓	✓	13
ENZO <sup>2</sup> 7 EN798	Super Power BTE	✓	✓	675
<b>ENZO<sup>2</sup> 5 Wireless</b>				
ENZO <sup>2</sup> 5 EN588	Power BTE	✓	✓	13
ENZO <sup>2</sup> 5 EN598	Super Power BTE	✓	✓	675

## CONSUMABLES

Metal Hook

ENZO<sup>2</sup> BTE's use Thick Tubing with Earmould.


## SOFTWARE AND PROGRAMMING

Aventa 3.10 or higher


Please see Cable and Wireless fitting guide on page 45 for details

## FREQUENTLY ASKED QUESTIONS


### Attaching the programming adaptor for ENZO<sup>2</sup> 98



Use a pin tool and tweezers to remove the battery door pin and the battery door.



Attach the programming adapter/door by inserting the pin again.



Insert a fresh battery and attach the CS44 cable.

# ReSound Up Smart

## DEVICE SUMMARY

MODEL	STYLE	MFI	WIRELESS ACCESSORIES COMPATIBILITY	BATTERY
<b>Up Smart 9 Wireless</b>				
Up Smart 9 UPS977	Standard BTE	✓	✓	13
Up Smart 9 UPS988	Power BTE	✓	✓	13
Up Smart 9 UPS998	Super Power BTE	✓	✓	675
<b>Up Smart 7 Wireless</b>				
Up Smart 7 UPS777	Standard BTE	✓	✓	13
Up Smart 7 UPS788	Power BTE	✓	✓	13
Up Smart 7 UPS798	Super Power BTE	✓	✓	675
<b>Up Smart 5 Wireless</b>				
Up Smart 5 UPS577M	Standard BTE	✓	✓	13
Up Smart 5 UPS588M	Power BTE	✓	✓	13
Up Smart 5 UPS598M	Super Power BTE	✓	✓	675

## CONSUMABLES

### Domes

Up Smart 77 Model



### Thin Tubes

Up Smart 77 Model



## SOFTWARE AND PROGRAMMING

Aventa 3.10 or higher

Please see Cable and Wireless fitting guide on page 45 for details

	Remote Control (1st Gen)	Phone Clip (1st Gen)	TV Streamer (1st Gen)	Mini Mic (1st Gen)	Remote Control 2 (2nd Gen)	Phone Clip+ (2nd Gen)	TV Streamer 2 (2nd Gen)	Micro Mic (2nd Gen)	Multi Mic (2nd Gen)
<b>GN ReSound</b>									
LiNX <sup>2</sup> 5   7   9	✓	✓	X	✓	✓	✓	✓	✓	✓
Enzo <sup>2</sup> 5   7   9	✓	✓	X	✓	✓	✓	✓	✓	✓
LiNX 5   7   9	✓	✓	✓	✓	✓*	✓	✓	✓	✓
Up Smart 5   7   9	✓	✓	✓	✓	✓	✓	✓	✓	✓
Enzo 5   7   9	✓	✓	✓**	✓	✓	✓	✓	✓	✓
Verso 5   7   9	✓	X	X	X	✓	X	X	X	X
Alera 5   7   9	✓	X	X	X	✓*	X	X	X	X
Enya 3   4	✓	✓	✓	✓	✓	✓	✓	✓	✓
Enya 2	✓	X	X	X	✓	X	X	X	X
Alera 4	✓	X	X	X	✓*	X	X	X	X
<b>Cochlear</b>									
BAHA 4	✓	✓	X	✓	✓	✓	✓	✓	✓
Nucleus 6	X	✓*	✓	✓	X	✓*	✓	✓	✓

**When mixing a Cochlear Implant and a Hearing Instrument, the Hearing Instrument must be paired first and the Implant is paired next, as soon as the pairing light starts to blink again.**

- \* Control part (program change, volume adjustment and mute functionality) not supported Alera/Nucleus.
- \*\* Not compatible with custom devices.